

SZ61/SZ51

Zoom Stereo Microscope

These compact stereo microscopes help you remain comfortable with eyepieces that reduce eyestrain, a universal LED stand that provides easy access to your sample, and a Greenough optical system that delivers excellent flatness. High color fidelity and antistatic materials and coatings make these microscopes suitable for routine or advanced life science applications.




SZ61/SZ51



SZ61TR



SZ61-60/SZ51-60



Observation Images

Specifications

	SZ61	SZ61-60	SZ61	SZ51	SZ51-60
Optical system	Greenough type optical system				
Zoom ratio	1:6.7			1:5	
Working distance	110 mm				
Tube inclination angle	45°	60°	45°		60°
Video camera adaptability	—		C-mount (0.5X built in)	—	
Optical component	Lead-free				
Auxiliary objective	Mounting by screwing into the thread at the bottom of frame (M48 threadx0.75)				
Eyepieces	WHSZ series				
Dimensions (W × D × H)	194 mm × 253 mm × 368 mm (Standard Set Configuration)				

SZX/SZ Accessories



SZ2-ILST

LED illuminator stand

This LED stand features a thin design to keep sample positions low and optimize operability. Simultaneous transmitted and reflected light are available on this stand. LED light offers a long lifespan and constant color temperature at any intensity.



SZ2-ILA

Transmitted illumination attachment

Used with the SZ2-ST stand, this cost-effective illumination stand offers bright, uniform illumination from low to high magnifications. A tiltable mirror provides direct and oblique illumination for low-contrast specimen. The available 6 W and 100 W lamphouses deliver the necessary power for various illumination needs.



SZ2-CLGDI

Dual interlock light guide

Standard oblique semi-rigid fiber optic light guide. The light source position on the rear side of the stand saves desk space.



SZ2-CLGSF

Flexible light guide

A single fiber optic guide is fixed at the back of the objective so as not to disturb microscope operation.



SZX2-ILLTQ/SZX2-ILLTS

LED Transmitted Light Illumination bases

SZX2-ILLTQ with quad-position turret enables you to choose cartridges and to switch between brightfield (standard/high/low), oblique (standard/high/low), darkfield, polarized illumination, and the shutter with a simple turn. A one-position LED illumination base is also an option (SZX2-ILLTS).



SZX-SDO2

Side-by-Side Discussion Tube

Ample distance (650 mm) between the main and secondary observer for easy imaging without disturbing the microscope operation. The color of the built-in pointer can be selected to contrast the specimen.

EVIDENT
OLYMPUS

LABEX
FLUOCCA
IMARIS

GREINER
BIO SCIENCE

APPLIED
SPECTRAL
IMAGING

SAKURA
FINETEK

ZTEC
CRAFTTEK

TEKSQRAY
FEATHER

DYNAMICA
FROLABO

FASTER
ARCTIKO
KÜHNER
HETTICH

HBK
CAMERA
IX CAMERAS

THERMO
EM
CRYSTAL