

EVIDENT
OLYMPUS

LABEX
FLUOCCA
IMARIS

GREINER
BIO SCIENCE

APPLIED
SPECTRAL
IMAGING

SAKURA
FINETEK

ZTEC
CRAFTTEK

TEKSQRAY
FEATHER

DYNAMICA
FROLILABO

FASTER
ARCHIKO
KUNHNER

IX CAMERAS






HBK
CAMECA
THERMO EM
CRYSTAL

SLIDEVIEW VS200

Research Slide Scanner

The SLIDEVIEW VS200 research slide scanner is designed to capture high-resolution images of your slides for quantitative analysis. The system offers even illumination and uses X Line high-performance objectives for better resolution and flatness. Fast batch scanning, a simplified workflow, and flexible imaging capabilities help you achieve more in less time.

Observation Images



Human cartilage captured with X Line UPLXAPO10X objective.

Specifications

Intended specimen	Observable Specimen	Glass slide with cover glass
	Size Options for Glass Slides	Width: 25 mm – 26 mm, length: 75 mm – 76 mm, thickness: 0.8 mm – 1.4 mm
	Size of Cover Glass	Thickness: 0.12 mm – 0.17 mm
Microscope frame	Illuminator	Built-in Kohler illumination for transmitted light
	Objectives (*option)	2X, 10X, 20X, 40X, 60X*1 and 100X*1 oil immersion with a motorized revolving nosepiece
	Motorized Stage	Motorized XY stage with automatic control
	Focusing	Motorized automatic control
Digital camera	Fluorescence Observation (Option)	Motorized fluorescence illuminator, motorized filter wheel, fluorescence light source, digital monochrome camera
	CCD Camera	2/3" CCD camera, 3.45 μ m $\times$ 3.45 μ m pixel size, high sensitivity, high resolution
	Image Correction	Shading correction, auto white balance
Scan	Scan Area	26 mm (W) $\times$ 64 mm (H) (Slide glass size: 26 mm (W) $\times$ 76 mm (H))
	Pixel Resolution *This value is length per pixel of the specimen surface.	20X (NA 0.75): 0.33 μ m/pixel 40X (NA 0.95): 0.17 μ m/pixel 60X oil immersion (option, NA 1.35): 0.11 μ m/pixel 100X oil immersion (option, NA 1.4): 0.07 μ m/pixel
	Scan Time	Approx. 2 min. (20X objective, scan area 15 mm $\times$ 15 mm)

