


		IX73	IX83
Microscope frame	Optical system	UIS2 optical system	
	Revolving nosepiece	- Motorized sextuple revolving nosepiece (DIC slider attachable), simple waterproof structure - Coded sextuple revolving nosepiece (DIC slider attachable), simple waterproof structure	Motorized sextuple revolving nosepiece (DIC slider attachable), simple waterproof structure
	Focus	Stroke: 10 mm	Stroke: 10.5 mm Minimum increment: 0.01 μm, Maximum nosepiece movement speed: 3 mm/s
	Intermediate port	2 ports (IX73P2F), 1 port (IX73PIF)	1 port (IX83PIZF)
	Light path selection	Manual 0.100/50.50/100.0 (Left side port: BI port)	Motorized 0.100/50.50/100.0 (Left side port: BI port)
Transmitted light illuminator		Pillar tilt mechanism (30° inclination angle, with vibration reducing mechanism), Condenser holder (with with 88 mm stroke, refocusing mechanism), Field iris diaphragm adjustable, 4 filter holders Light source: 12 V, 100 W halogen bulb (pre-centered) - High color reproductive LED light source	
Observation tube	Widefield (FN 22)	- Widefield tilting binocular - Widefield binocular	
Stage	Mechanical stage with left handle	Stage stroke: X: 50 mm x Y: 43 mm, including stage insert plate (ø25 mm), stage insert plate exchangeable (Ø110 mm)	
	Cross stage with right handle	Stage stroke: X: 50 mm x Y: 50 mm, including stage insert plate (ø25 mm), stage insert plate exchangeable (Ø110 mm)	
	Mechanical stage with right handle	Stage stroke: X: 114 mm x Y: 75 mm, stage position locking function	
Condenser	Motorized middle long working distance condenser	NA 0.5, W.D. 45 mm, 4 positions for optical devices (for ø50 mm, relief contrast optical devices rotatable)	
	Mid long working distance relief contrast	NA 0.5, W.D. 45 mm, 4 positions for optical devices (for ø50 mm, relief contrast optical devices rotatable)	
Operating environment	- Indoor use - Ambient temperature: 5° to 40°C (41° to 104°F) - Maximum relative humidity: 80% for temperatures up to 31°C (88°F), decreasing linearly through 70% at 34°C (93°F), 60% at 37°C (99°F), to 50% relative humidity at 40°C (104°F) - Supply voltage fluctuations: Not to exceed ±10% of the normal voltage		

Morphological Evaluation of Sperm at High Magnification

Improve the sperm selection process:

- Smooth zooming from 1x, 1.6x, to 2x using the adjustment lever.
- Combine the zoom with the 40X objective lens to observe the morphology and internal structure of the sperm head in detail.



Learn More

