

IX Applications for ICSI/IMSI

Optimized for a Smoother ICSI Workflow

IX73 Microscope

Motorized Observation Method Switching

- Quickly change observation methods by pressing a button on the hand switch
- Spindle observation is standard with a semi-motorized system; IMSI units are available




IX83 Microscope

Fully Motorized and Highly Advanced

- Standard configuration supports IMSI and spindle observation
- Quickly change observation methods and magnification using the controller
- Focus position and brightness automatically adjust when you change magnifications
- Easily operate the focus and switch between camera and visual observation using the remote control, reducing sample vibration



Applied Technologies




Easy Oocyte Condition Checks through Spindle Visualization

Olympus' observation method enables you to easily view the metaphase II (MII) spindle through the microscope oculars



Relief Contrast

Olympus' relief contrast enables 3D imaging of oocytes so you can check the condition of the zona pellucida



Polarization

Observe the birefringence of spindles using polarized light